# **VEVOR**


**Upgrade · The Home Creator Way**

## **LAB CHILLER CIRCULATOR**

**MODEL:DC-0506/DC-2006/DC-3006**

**MODEL:DC-0506/DC-2006/DC-3006**

This is the original instruction, please read all manual instructions carefully before operating. VEVOR reserves a clear interpretation of our user manual. The appearance of the product shall be subject to the product you received. Please forgive us that we won't inform you again if there are any technology or software updates on our product.

	<p>Warning-To reduce the risk of injury, the user must read the instructions manual carefully.</p>
	<p>This product is subject to the provision of European Directive 2012/19/EC. The symbol showing a wheeled bin crossed through indicates that the product requires separate refuse collection in the European Union. This applies to the product and all accessories marked with this symbol. Products marked as such may not be discarded with normal domestic waste, but must be taken to a collection point for recycling electrical and electronic devices</p>


## SAFETY INSTRUCTIONS AND PRECAUTIONS

### WARNING:

1. Before using, the slot should be in a liquid medium. (Note: liquid does not contain acid or alkali).
2. The total power of the power is to be greater than or equal to the instrument. The power must have a good "grounding" device. Using the voltage will be subject to technical parameter table.
3. The equipment should be placed in a well ventilated and dry place; Do not leave any obstacles on the back and side within a distance of 300mm.
4. After use, switch off and unplug the power plug.

# Instrument appearance introduction

1	Control box	2	The instruent set key
3	The shift key	4	The reduction key
5	The addend key	6	The circular key
7	Refrigeration	8	Power key
9	Power Switch	10	Power Socket
11	Water Inlet	12	Water Outlet
13	Revolving Insurance		

## Features

1 Refrigeration System: Low-noise air-cooled totally enclosed compressor unit is adopted, which has the advantages of fast refrigeration speed and strong temperature stability.

2 Control System: Latest proprietary software research and development of homemade temperature control technology, with high standards of PT100 and all imports of electronic parts, etc.

3 Circulating System: It is composed of internal and external circulation. During external circulation, the constant temperature liquid in the tank can be introduced outside, the second constant temperature field can be established, and it can also be used as a cold source to lead the cooled liquid in the tank to the experimental

container outside the machine for cooling; During internal circulation, the liquid temperature in the tank can be very uniform and stable.

4 Protection System: Refrigeration is overheated, over-current protection, call control system has over temperature alarm, can be set over temperature alarm temperature, the temperature of the upper and lower limits can cut off the load automatically when over temperature, and other functions.

5 High-end LCD control instrument, simple operation, temperature stability is strong, P.i.d can automatically search or artificial adjustment

6 The accurate temperature correction function can reach  $0.01^{\circ}\text{C}$ , and the maximum temperature fluctuation can reach  $\pm 0.02-0.05^{\circ}\text{C}$  (depending on the model).


## Technique Parameter

Model	DC-0506	DC-2006	DC-3006
Temperature range ( $^{\circ}\text{C}$ )	-5~100	-20~100	-30~100
Temperature fluctuation ( $\pm^{\circ}\text{C}$ )	0.03	0.03	0.03
Groove depth(mm)	150	150	150
Power (W)	120V~ 60Hz 1000W	120V~ 60Hz 1000W	120V~ 60Hz 950W
	230V~ 50Hz 800W	230V~ 50Hz 900W	230V~50Hz 1050W

## Operating Steps


1. Connect the circulating pump before startup (if there is no water inlet and outlet, this step can be omitted):

1.1 For the connection of circulating pump during internal circulation, connect the water outlet and water inlet directly with both ends of the hose supplied with the goods.

Simply connect the water outlet and inlet with both ends of the hose provided in the package.


1.2 For the connection of circulating pump during external circulation, connect the water outlet to the inlet of container / equipment outside the tank with hose, and connect the water inlet to the outlet of container / equipment outside the tank.

Connect the water outlet to the inlet of the container/ equipment outside the tank with a hose.

Connect the water inlet to the outlet of the container/ equipment outside the tank.

3. Add the corresponding liquid medium into the tank. The liquid level of the liquid medium shall not be 20mm lower than the worktable to avoid damage caused by exposed dry burning of the heater.

2.1 When the operating temperature is below 8°C, industrial alcohol or antifreeze is generally used for liquid media.

2.2 At an operating temperature of 8°C-75°C, pure water is generally selected for liquid media.

2.3 At an operating temperature of 75°C-100°C, dimethyl silicone oil with a viscosity of 5cs is generally selected for liquid media.

3.Drain the liquid in the tank.


When the equipment is not used for a long time or the liquid is replaced, the liquid in the tank needs to be drained.

3.1Do not operate under high-temperature conditions to avoid scalding.

3.2When the equipment is out of use for an extended period or when replacing the liquid, the liquid in the tank must be drained.



3.3When the equipment is shut down and the ambient temperature is below 4°C, the water in the unit's tank must be completely drained to prevent freezing and cracking of the water tank and pipelines.

3.4Connect the drain nipple with a hose, then open the valve to let the liquid drain out.


Attach a hose to the liquid drain barb fitting, then activate the valve to allow liquid evacuation.

4. Plug in the power, turn on the "main power" switch, and then turn on the "power" button on the operation panel. (Note: "main power" switch is on the back of the instrument).

5. Perform the instrument operation as follows:

5.1 Display Description:


Icon	Name	Button Description
	Set Key	Click to set the temperature or long - press this key to enter the temperature parameter setting state.
	Shift Key	In the setting state, click this key to make the set value shift and blink for modification; in the main interface display state, long - press this key for 6 seconds to enter the temperature self - calibration selection and self - calibration value setting state.
	Decrease Key	In the setting state, click or long - press this key to decrease the set value.
	Increase Key	In the setting state, click or long - press this key to increase the set value.
	Power Key	Turn off or turn on the controller.
	Cooling Key	Click to allow or prohibit cooling
	Cycle Key	Click to turn on or turn off the cycle output.
	Overtemperature Alarm	The difference between the set temperature and the actual temperature exceeds 30°C
	Timer Function	0000 means no timer is set and the device works all the time; 0001 - 9999min means the operable total running time after the device reaches the set temperature point. When the time is up, the device stops running and the character "END" appears
E-01	Alarm	Phase Sequence Alarm
E-02	Alarm	Liquid Level Alarm

## 5.2 Description of instrument keys:

Set function key shift key ▼ decrease key ▲ increase key

## 5.3 Setting of temperature parameters:

Press the setting function key "set" to enter the setting state of the temperature setting value. At this time, the last value of the setting temperature "SV" on the display screen flashes.

First, press the shift key, then press the plus or minus key to set the working temperature you need, and then press the setting function key "set" once. At this time, the last value of the timed shutdown time (min) "0000" on the upper right corner of the display screen flashes. If it is not necessary to keep the default "0000" regularly, press the setting function key "set" again to exit the parameter setting and save the parameter setting value(PV). At this time, the display "PV" shows the temperature of the liquid medium in the current tank.

5.4 After the temperature parameter setting is completed, first open the "circulate" button on the operation panel, and then open the "refrigeration" button (Note: when the working temperature exceeds 40 °C, the compressor automatically enters the protection, and the refrigeration does not start), and then the microcomputer enters the automatic control state.

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