

PROCAPMED



STANDS ON DURABILITY,
WORKS ON RELIABILITY

Procapmed Pipettes

Micro Volume Pipette

Product Features



Ergonomic &
Light Weight Design



Smooth Plunger Mechanism
Offers Least Resistance For
Ease-of-use And Long Hours
Of Pipetting



High Accuracy And
Precision Guaranteed.




Durable Tip Cone Made Of
PVDF (Polyvinylidene Fluoride)
Offers The Highest Resistance
To Physical And Chemical
Corrosion And Is Compatible
With Wide Range Of Tips

Tip Ejector For Smooth
One-hand Ejection Of
The Tip Preventing
Any Contamination.



Two Step Plunger
Operation Allows Forward &
Reverse Pipetting Technique



Easy Click Stop
Volume Setting



Procapmed Pipettes

Micro Volume Pipette

Precision and accuracy are of the utmost importance when it comes to using pipettes, but existing products tend to fall short as they don't deliver reliable precision. In addition, most of these pipettes are not fully autoclavable which leads to greater inconvenience for users.


After years spent researching about customers' expectations and requirements related to micropipettes, Procapmed pipettes which combine a fully autoclavable setup with a lightweight and ergonomic body for a hassle and pain-free pipetting experience while maintaining unmatched accuracy and precision.



Product Features



Ideal Shaft Shape Helps
The Pipette Reach Contents
In Narrow Or Deep Tubes.



Fully Autoclavable



Calibrated As Per
ISO 8655 In ISO 17025
Accredited Lab

3 Year
Extended Warranty



Individual Serial No.
Permanently Marked
On Each Pipette For
Easy Traceability



Multichannel Pipettes

Includes All Basics Features Of Procapmed Single Channel & Other Feature Like

- Soft force / leak free tip loading and ejection.
- Easily removable & fully autoclavable manifold.
- Flexibility to operate pipette without twisting the wrist with the help of 360° rotatable Manifold.
- One hand tip ejection
- Attractive color coding
- Compatible with wide range of tips
- Available in 8 & 12 channels



Ordering Information

Procapmed Fully Autoclavable Variable Volume Micropipette

Cat. No.	Volume Range [µl]	Increments [µl]	Test Volume [µl]	Inaccuracy (±) %	Inaccuracy (±) µl	Imprecision (±) %	Imprecision (±) µl
CV-100	0.5-10	0.1	1	2.5	0.03	1.5	0.02
			5	2	0.10	1	0.05
			10	1	0.10	0.8	0.08
CV-125	0.1 - 2.5	0.01	0.25	12	0.03	6	0.02
			1.25	3	0.04	3	0.04
			2.5	2.5	0.06	1.6	0.04
CV-200	5-50	0.5	5	2	0.10	2	0.10
			25	0.8	0.20	0.4	0.10
			50	0.6	0.30	0.3	0.15
CV-500	10-100	0.5	10	3	0.30	1.5	0.15
			50	1	0.50	0.5	0.25
			100	0.8	0.80	0.15	0.15
CV-600	100-1000	5	100	2	2.00	0.7	0.70
			500	1	5.00	0.4	2.00
			1000	0.6	6.00	0.2	2.00
CV-700	2-20	0.1	2	3	0.06	2	0.04
			10	1.5	0.15	1	0.10
			20	0.9	0.18	0.4	0.08
CV-900	1000-5000	50	1000	0.8	8.00	0.3	3.00
			2500	0.6	15.00	0.2	5.00
			5000	0.5	25.00	0.15	7.50
CV-1000	1000-10000	100	1000	3	30.00	0.6	6.00
			5000	1.2	60.00	0.3	15.00
			10000	0.6	60.00	0.2	20.00

Ordering Information

Procapped Fully Autoclavable Fix Volume Micropipette						
Cat. No.	Volume [μ l]	Test Volume	Inaccuracy ($\pm$) %	Inaccuracy ($\pm$) μ l	Imprecision ($\pm$) %	Imprecision ($\pm$) μ l
CF-2.5	2.5	2.5	2	0.05	1.6	0.04
CF-5	5	5	1.3	0.07	1.2	0.06
CF-10	10	10	0.8	0.08	0.8	0.08
CF-20	20	20	0.6	0.12	0.5	0.10
CF-25	25	25	0.5	0.13	0.3	0.08
CF-50	50	50	0.5	0.25	0.3	0.15
CF-100	100	100	0.5	0.50	0.3	0.30
CF-200	200	200	0.4	0.80	0.2	0.40
CF-250	250	250	0.4	1.00	0.2	0.50
CF-500	500	500	0.3	1.50	0.2	1.00
CF-1000	1000	1000	0.3	3.00	0.2	2.00
CF-2000	2000	2000	0.3	6.00	0.15	3.00

Ordering Information

Multi Channel Pipettes - 8							
Cat. No.	Volume [µl]	Increments [µl]	Test Volume [µl]	Inaccuracy (±) %	Inaccuracy (±) µl	Imprecision (±) %	Imprecision (±)µl
CM8-010	0.5-10	0.1	10	1.5	0.15	1.5	0.15
			5	2.5	0.13	2.5	0.13
			1	4	0.04	4	0.04
CM8-050	5-50	0.5	50	1	0.50	0.7	0.35
			25	1.5	0.38	1	0.25
			5	3	0.15	2	0.10
CM8-100	10 – 100	1	100	1	1.00	0.5	0.50
			50	1	0.50	0.5	0.25
			10	1.5	0.15	0.75	0.08
CM8-200	20-200	1	200	0.7	1.40	0.25	0.50
			100	1	1.00	0.4	0.40
			20	1.5	0.30	0.75	0.15
CM8-300	30-300	1	300	0.8	2.40	0.25	0.75
			150	1	1.50	0.5	0.75
			30	1.5	0.45	0.75	0.23